


Anti-Tcp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17987

Item Number: CM0015900



Introduction

Rabbit polyclonal antibody targeting human TCP1 (CCT-alpha), validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat TCP1. Ideal for chaperonin, cytoskeleton, and ciliogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	TCP1 Rabbit pAb
Remarks/Alias	CCT1; CCTa; D6S230E; CCT-alpha; TCP-1-alpha; TCP1
Species	Human
GeneID (Human)	6950
GeneID	6950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-556 of human TCP1 (NP_110379.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17987
Protein Weight	60kDa
Shipping	Ice bag