

Anti-Tcta Polyclonal Antibody For Wb, If/Icc, Elisa - P57738

Item Number: CM0014143

Introduction

TCTA Rabbit pAb (SWISS: P57738) is a polyclonal antibody raised against full-length human TCTA protein. Validated for WB, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for osteoclast and membrane protein studies.

[Learn More](#)

Parameter	Value
Product Name	TCTA Rabbit pAb
Remarks/Alias	TCTA
Species	Human
GeneID (human)	6988
GeneID	6988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human TCTA (NP_071503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57738
Protein Weight	11kDa
Shipping	Ice bag