


Anti-Tcte1 Polyclonal Antibody For Wb, Elisa - Q5Ju00

Item Number: CM0015256



Introduction

Rabbit polyclonal antibody against human TCTE1 (DRC5), suitable for WB and ELISA. Recognizes human, mouse, and rat TCTE1. Useful for investigating ciliary/flagellar motility and N-DRC function.

[Learn More](#)

Parameter	Value
Product Name	TCTE1 Rabbit pAb
Remarks/Alias	DRC5; D6S46; FAP155; TCTE1
Species	Human
GeneID (Human)	202500
GeneID	202500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TCTE1 (NP_872345.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5JU00
Protein Weight	56kDa
Shipping	Ice bag