


Anti-Tctp/Tpt1 Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P13693

Item Number: CM0006830



Introduction

High-quality Anti-TCTP/TPT1 Rabbit Monoclonal Antibody validated for WB, IP, IHC-P, ELISA in Human, Mouse, Rat. Recognizes 20kDa protein involved in calcium binding, microtubule stabilization, and apoptosis regulation. IVD raw material for research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	TCTP/TPT1 Rabbit mAb
Remarks/Alias	HRF; p02; p23; TCTP; TCTP/TPT1
Species	Human
GeneID (human)	7178
GeneID	7178
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TCTP/TPT1 (P13693)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13693
Protein Weight	20kDa
Shipping	Ice bag