


Anti-Tctp/Tpt1 Rabbit Polyclonal Antibody For Wb And Elisa - P13693

Item Number: CM0027395



Introduction

Rabbit polyclonal antibody targeting human TCTP/TPT1, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for calcium signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TCTP/TPT1 Rabbit pAb
Remarks/Alias	HRF; p02; p23; TCTP; TCTP/TPT1
Species	Human
GeneID (Human)	7178
GeneID	7178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human TCTP/TCTP/TPT1 (NP_003286.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13693
Protein Weight	20kDa
Shipping	Ice bag