


Anti-Tctp/Tpt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13693

Item Number: CM0013801



Introduction

Rabbit polyclonal antibody targeting human TCTP/TPT1 (HRF, p23, fortilin). Recombinant immunogen (aa 1-172). Suited for WB, IF/ICC, ELISA; cross-reacts with human, mouse. Supports research on apoptosis, calcium binding, and microtubule dynamics.

[Learn More](#)

Parameter	Value
Product Name	TCTP/TPT1 Rabbit pAb
Remarks/Alias	HRF, p02, p23, TCTP, TCTP/TPT1
Species	Human
GeneID (Human)	7178
GeneID	7178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human TCTP/TCTP/TPT1 (NP_003286.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13693
Protein Weight	20kDa
Shipping	Ice bag