

Anti-Tdg Polyclonal Antibody For Wb And Elisa - Q13569

Item Number: CM0026129



Introduction

Rabbit polyclonal antibody against human TDG, validated for Western blot and ELISA. Suitable for human, mouse, and rat samples. Ideal for studies of DNA demethylation and base excision repair.

[Learn More](#)

Parameter	Value
Product Name	TDG Rabbit pAb
Remarks/Alias	hTDG; TDG
Species	Human
GeneID (Human)	6996
GeneID	6996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TDG (NP_003202.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13569
Protein Weight	46kDa
Shipping	Ice bag