

Anti-Tdgf1/Cripto Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P13385

Item Number: CM0018383

Introduction

Rabbit polyclonal antibody against human TDGF1 (CRIPTO), validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat TDGF1. Suitable for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	TDGF1 Rabbit pAb
Remarks/Alias	CR; CR-1; CRGF; CRIPTO; TDGF1
Species	Human
GenelD (Human)	6997
GenelD	6997
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-150 of human TDGF1 (NP_001167607.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13385
Protein Weight	21kDa
Shipping	Ice bag