


Anti-Tdo2 Polyclonal Antibody For Wb, If/Icc, Elisa - P48775

Item Number: CM0031041



Introduction

This anti-TDO2 rabbit polyclonal antibody detects human and mouse tryptophan 2,3-dioxygenase (~48 kDa) and is validated for WB, IF/ICC, and ELISA applications. Ideal for studies on tryptophan catabolism and related diseases.

[Learn More](#)

Parameter	Value
Product Name	TDO2 Rabbit pAb
Remarks/Alias	TO; TDO; TPH2; TRPO; HYPTRP; TDO2
Species	Human
Gene ID (Human)	6999
Gene ID	6999
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TDO2 (NP_005642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48775
Protein Weight	48kDa
Shipping	Ice bag