

Anti-Tdp-43/Tardbp Monoclonal Antibody For Wb, Ip, If, Ihc, Chip - Q13148

Item Number: CM0003765



Introduction

Recombinant rabbit monoclonal antibody against human TDP-43/TARDBP (Q13148). Suitable for WB, IP, IF, IHC, ChIP, ELISA. Cross-reacts with mouse and rat. Evaluated for research on RNA processing and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	TDP-43/TARDBP Rabbit mAb
Remarks/Alias	ALS10; TDP-43; TDP-43/TARDBP
Species	Human
GeneID (human)	23435
GeneID	23435
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 300-400 of human TDP-43/TARDB (Q13148))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13148
Protein Weight	45 kDa
Shipping	Ice bag