

Anti-Tdp-43/Tardbp Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13148

Item Number: CM0017958



Introduction

Rabbit polyclonal antibody against human TDP-43/TARDBP, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for RNA biology, neurodegeneration, and circadian rhythm research.

[Learn More](#)

Parameter	Specification
Product Name	TDP-43/TARDBP Rabbit pAb
Remarks/Alias	ALS10; TDP-43; TDP-43/TARDBP
Species	Human
GeneID (human)	23435
GeneID	23435
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13148
Protein Weight	45kDa
Shipping	Ice bag