

Anti-Tdp-43/Tardbp Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa, Rip - Q13148

Item Number: CM0017850



Introduction

Rabbit polyclonal antibody against human TDP-43/TARDBP (Q13148), validated for WB, IHC-P, IP, ELISA, and RIP. Cross-reacts with mouse and rat. Suitable for neurodegenerative disease and RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	TDP-43/TARDBP Rabbit pAb
Remarks/Alias	ALS10; TDP-43; TDP-43/TARDBP
Species	Human
GeneID (Human)	23435
GeneID	23435
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TDP-43/TARDBP (NP_031401.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA, RIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13148
Protein Weight	45kDa
Shipping	Ice bag