

Anti-Tdp1 Rabbit Pab For Wb, If/Icc, Elisa - Q9Nuw8

Item Number: CM0029956



Introduction

Rabbit polyclonal antibody recognizing human and rat TDP1 (68 kDa). Validated for WB, IF/ICC, ELISA; immunogen covers amino acids 1-140. Suited for DNA repair and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	TDP1 Rabbit pAb
Remarks/Alias	TDP1
Species	Human
Gene ID (Human)	55775
Gene ID	55775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human TDP1 (NP_060789.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NUW8
Protein Weight	68kDa
Shipping	Ice bag