


Anti-Tdp2 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - O95551

Item Number: CM0020264



Introduction

Rabbit polyclonal antibody against TDP2, KO validated. Detects human and mouse TDP2 (O95551) via WB and ELISA. Suitable for DNA repair, topoisomerase biology, and neurological research applications.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TDP2 Rabbit pAb
Remarks/Alias	EAP2; AD022; EAPII; TTRAP; hTDP2; dj30M3.3; P2
Species	Human
GeneID	51567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TDP2 (NP_057698.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95551
Protein Weight	41kDa
Shipping	Ice bag