

Anti-Tdrd1 Polyclonal Antibody For Wb, Elisa - Q9Bxt4

Item Number: CM0015340



Introduction

Rabbit polyclonal TDRD1 antibody validated for Western blot and ELISA. Cross-reacts with mouse. Targets the cancer/testis antigen CT41.1, essential for piRNA-mediated transposon silencing in spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	TDRD1 Rabbit pAb
Remarks/Alias	CT41.1; TDRD1
Species	Human
GeneID	56165
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 25-125 of human TDRD1 (NP_942090.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BXT4
Protein Weight	132kDa
Shipping	Ice bag