

Anti-Tdrd3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H7E2

Item Number: CM0025932



Introduction

Rabbit polyclonal antibody against human TDRD3, validated for WB and ELISA. Recognizes a scaffolding protein involved in transcriptional coactivation and stress granule assembly. Research use only.

[Learn More](#)

Parameter	Value
Product Name	TDRD3 Rabbit pAb
Remarks/Alias	TDRD3
Species	Human
Gene ID (Human)	81550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TDRD3 (NP_110421.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H7E2
Protein Weight	73kDa
Shipping	Ice bag