


Anti-Tdrkh Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2W6

Item Number: CM0027650



Introduction

High-quality rabbit polyclonal antibody against TDRKH, suitable for WB and ELISA. Shows cross-reactivity with human, mouse, and rat. Ideal for investigating piRNA biogenesis and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	TDRKH Rabbit pAb
Remarks/Alias	TDRD2; TDRKH
Species	Human
GeneID (human)	11022
GeneID	11022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-355 of human TDRKH (NP_001077432.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2W6
Protein Weight	62kDa
Shipping	Ice bag