


Anti-Tdt Rabbit Polyclonal Antibody For Wb, Elisa - P04053

Item Number: CM0029498



Introduction

Rabbit polyclonal antibody targeting human TdT, validated for WB and ELISA. Reacts with human, mouse, rat. Useful for studying V(D)J recombination and immune cell development.

[Learn More](#)

Parameter	Value
Product Name	TdT Rabbit pAb
Remarks/Alias	TDT; TdT
Species	Human
GeneID (Human)	1791
GeneID	1791
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04053
Protein Weight	59kDa
Shipping	Ice bag