


Anti-Tead1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - P28347

Item Number: CM0026661



Introduction

Rabbit polyclonal antibody against TEAD1, validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Cross-reacts with human, mouse, rat. Ideal for Hippo signaling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	TEAD1 Rabbit pAb
Remarks/Alias	AA; REF1; TCF13; TEF-1; NTEF-1; TCF-13; TEAD-1; TEAD1
Species	Human
Gene ID (Human)	7003
Gene ID	7003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-215 of human TEAD1 (NP_068780.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P28347
Protein Weight	48kDa
Shipping	Ice bag