

Anti-Tead1/2/3/4 Rabbit Monoclonal Antibody For Wb, Elisa - P28347 / Q15561 / Q15562 / Q99594

Item Number: CM0012743



Introduction

High-quality anti-TEAD1/2/3/4 rabbit monoclonal antibody for Western blotting and ELISA, targeting human TEAD transcription factors in Hippo signaling research.

[Learn More](#)

Parameter	Specification
Product Name	TEAD1/2/3/4 Rabbit mAb
Remarks/Alias	AA; REF1; TCF13; TEF-1; NTEF-1; TCF-13; TEAD-1
Species	Human
Gene ID (Human)	7003/7004/8463/7005
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 174-263 of human TEAD4 (NP_003204.2)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P28347/Q15561/Q15562/Q99594
Protein Weight	48 kDa/49 kDa/47 kDa
Shipping	Ice bag