

Anti-Tead2 Rabbit Polyclonal Antibody For Wb, Elisa - Q15562

Item Number: CM0012961



Introduction

Rabbit polyclonal antibody targeting human TEAD2 (TEF-4), validated for WB and ELISA. Detects ~49 kDa protein. Essential for Hippo pathway, cell proliferation, and transcriptional regulation studies.

[Learn More](#)

Specification	Detail
Product Name	TEAD2 Rabbit pAb
Remarks/Alias	ETF; TEF4; TEF-4; TEAD-2; TEAD2
Species	Human
GeneID (Human)	8463
GeneID	8463
Immunogen	Recombinant protein of human TEAD2
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15562
Protein Weight	49kDa
Shipping	Ice bag