

Anti-Tead4 Rabbit Monoclonal Antibody For Wb, Elisa - Q15561

Item Number: CM0006755



Introduction

Anti-TEAD4 rabbit monoclonal antibody for Western blot (WB) and ELISA to detect human TEAD4, a key Hippo signaling transcription factor regulating cell proliferation and migration.

[Learn More](#)

Parameter	Value
Product Name	TEAD4 Rabbit mAb
Remarks/Alias	TEAD4; EFTR-2; RTEF1; TCF13L1; TEF-3; TEF3; TEFR-1; hRTEF-1B; TEA domain transcription factor 4
Species	Human
Gene ID (Human)	7004
Gene ID	7004
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 140-240 of human TEAD4 (Q15561)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15561
Protein Weight	34kDa/44kDa/48kDa
Shipping	Ice bag