

Anti-Tead4 Rabbit Polyclonal Antibody For Wb, Elisa - Q15561

Item Number: CM0030880



Introduction

Rabbit polyclonal antibody targeting human TEAD4 (TEF-3), validated for WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for Hippo pathway research.

[Learn More](#)

Parameter	Value
Product Name	TEAD4 Rabbit pAb
Remarks/Alias	TEF3; RTEF1; TEF-3; EFTR-2; TEFR-1; TCF13L1; hRTEF-1B; TEAD4
Species	Human
Gene ID (Human)	7004
Gene ID	7004
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 110-250 of human TEAD4 (NP_003204.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15561
Protein Weight	48kDa
Shipping	Ice bag