


Anti-Tec Monoclonal Antibody For Wb And Elisa - P42680

Item Number: CM0004839



Introduction

High-quality rabbit monoclonal antibody targeting human Tec kinase. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for immune signaling research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Tec Rabbit mAb
Remarks/Alias	PSCTK4; Tec
Species	Human
GeneID	7006
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human Tec (NP_003206.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42680
Protein Weight	74kDa
Shipping	Ice bag