


Anti-Tec Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P42680

Item Number: CM0013148



Introduction

Tec Rabbit pAb is a polyclonal antibody raised against recombinant human Tec (aa 1-280). Validated in WB, IF/ICC, IHC-P, and ELISA, it detects Tec in human, mouse, and rat samples, making it ideal for immune signaling and kinase research.

[Learn More](#)

Parameter	Value
Product Name	Tec Rabbit pAb
Remarks/Alias	PSCTK4; Tec
Species	Human
GeneID (Human)	7006
GeneID	7006
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human Tec (NP_003206.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42680
Protein Weight	74 kDa
Shipping	Ice bag