


Anti-Tefm Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Qe5

Item Number: CM0028647



Introduction

High-quality rabbit polyclonal antibody against human TEFM (C17orf42), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying mitochondrial transcription and oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	TEFM Rabbit pAb
Remarks/Alias	C17orf42; TEFM
Species	Human
Gene ID (Human)	79736
Gene ID	79736
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-210 of human TEFM (NP_078959.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96QE5
Protein Weight	42kDa
Shipping	Ice bag