


Anti-Tek Rabbit Polyclonal Antibody For Wb And Elisa - Q02763

Item Number: CM0018027



Introduction

Rabbit polyclonal antibody against human TEK (TIE2, CD202B) for Western Blot and ELISA applications. Detects human, mouse, and rat TEK with a predicted band of 126 kDa. Suitable for angiogenesis and vascular biology studies.

[Learn More](#)

Parameter	Value
Product Name	TEK Rabbit pAb
Remarks/Alias	TIE2; VMCM; GLC3E; TIE-2; VMCM1; CD202B; TEK
Species	Human
GeneID (Human)	7010
GeneID	7010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-750 of human TEK (NP_000450.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02763
Protein Weight	126kDa
Shipping	Ice bag