


Anti-Tek (Tie2) Rabbit Polyclonal Antibody For Wb, Elisa - Q02763

Item Number: CM0024926



Introduction

High-quality rabbit polyclonal antibody targeting human TEK (TIE2; CD202B), validated for Western blot and ELISA. Suitable for angiopoietin receptor studies in angiogenesis and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	TEK Rabbit pAb
Remarks/Alias	TIE2; VMCM; GLC3E; TIE-2; VMCM1; CD202B; TEK
Species	Human
GeneID (Human)	7010
GeneID	7010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-750 of human TEK (NP_000450.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02763
Protein Weight	126kDa
Shipping	Ice bag