

Anti-Telo2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4R8

Item Number: CM0022394



Introduction

High-quality rabbit polyclonal antibody against human TELO2 (hCLK2, TEL2), validated for WB and ELISA. Reacts with human and mouse. Ideal for DNA damage response, mTOR signaling, and telomere biology research.

[Learn More](#)

Parameter	Value
Product Name	TELO2 Rabbit pAb
Remarks/Alias	CLK2; TEL2; YHFS; TELO2
Species	Human
Gene ID (Human)	9894
Gene ID	9894
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 450-550 of human TELO2 (NP_057195.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y4R8
Protein Weight	92kDa
Shipping	Ice bag