

Anti-Tenascin-X (Tnxb) Rabbit Polyclonal Antibody For Wb, Elisa - P22105

Item Number: CM0019985



Introduction

Rabbit polyclonal antibody targeting human Tenascin-X (TNXB) validated for WB and ELISA. Demonstrates cross-reactivity with mouse and rat. Suitable for extracellular matrix and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	TNXB Rabbit pAb
Remarks/Alias	XB; TNX; XBS; EDS3; HXBL; TENX; TN-X; VUR8; TNXB1; TNXB2; TNXB5; EDSCLL; EDSCLL1; TNXB
Species	Human
GeneID	7148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 444-673 of human TNXB (NP_115859.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22105
Protein Weight	458kDa
Shipping	Ice bag