

Anti-Tenm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukz4

Item Number: CM0015413

Introduction

Rabbit polyclonal antibody targeting human TENM1 (Teneurin-1). Validated for WB, IF/ICC, ELISA applications in human, mouse, and rat samples. Ideal for neural development and neuroplasticity research.

[Learn More](#)

Parameter	Value
Product Name	TENM1 Rabbit pAb
Remarks/Alias	TNM; ODZ1; ODZ3; TEN1; TNM1; ten-1; TEN-M1; TENM1
Species	Human
GeneID (human)	10178
GeneID	10178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-540 of human TENM1 (NP_001156750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKZ4
Protein Weight	305kDa
Shipping	Ice bag