


Anti-Terf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P54274

Item Number: CM0013162



Introduction

High-quality rabbit polyclonal antibody against human TERF1. Validated for WB, IHC-P, and ELISA. Cross-reacts with human and rat. Ideal for telomere biology, chromosome stability, and shelterin complex research.

[Learn More](#)

Parameter	Specification
Product Name	TERF1 Rabbit pAb
Remarks/Alias	TRF; PIN2; TRF1; TRBF1; t-TRF1; hTRF1-AS; TERF1
Species	Human
Gene ID	7013 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-275 of human TERF1 (NP_059523.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P54274
Protein Weight	50kDa
Shipping	Ice bag