


## Anti-Terf2 Polyclonal Antibody For Wb - Q15554

Item Number: CM0023114



### Introduction

Rabbit polyclonal antibody against human TERF2, validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for telomere biology and shelterin complex research.

[Learn More](#)

| Parameter        | Specification                                                                                         |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | TERF2 Rabbit pAb                                                                                      |
| Remarks/Alias    | TRF2; TRBF2; TERF2                                                                                    |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 7014                                                                                                  |
| GeneID           | 7014                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 200-300 of human TERF2 (NP_005643.2) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q15554                                                                                                |
| Protein Weight   | 60kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |