


Anti-Terf2Ip Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nyb0

Item Number: CM0026056



Introduction

TERF2IP (RAP1) rabbit polyclonal antibody for WB, IF/ICC, ELISA, validated in human and mouse. Targets telomeric repeat-binding factor 2-interacting protein 1, a shelterin component and NF- κ B regulator.

[Learn More](#)

Parameter	Value
Product Name	TERF2IP Rabbit pAb
Remarks/Alias	RAP1; DRIP5; TERF2IP
Species	Human
Gene ID	54386 (Human)
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TERF2IP (NP_061848.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NYB0
Protein Weight	44kDa
Shipping	Ice bag