


## Anti-Tert Monoclonal Antibody For Wb, Elisa - O14746

Item Number: CM0012515



### Introduction

High-quality rabbit monoclonal antibody targeting human TERT (Telomerase reverse transcriptase), validated for WB and ELISA. Cross-reacts with mouse and rat. Essential tool for telomere biology, aging, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TERT Rabbit mAb                                                                                                                       |
| Remarks/Alias    | TP2; TRT; CMM9; EST2; TCS1; hTERT; DKCA2; DKCB4; hEST2; PFBMFT1; TERT                                                                 |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 7015                                                                                                                                  |
| GeneID           | 7015                                                                                                                                  |
| Immunogen        | Recombinant protein   Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1132 of human TERT (O14746). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | O14746                                                                                                                                |
| Protein Weight   | 127kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |