

Anti-Tert Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14746

Item Number: CM0020901



Introduction

Rabbit polyclonal antibody targeting human TERT (telomerase reverse transcriptase). Validated for WB, IHC-P, ELISA. Reacts with human and mouse. Ideal for telomerase, cancer, and aging research.

[Learn More](#)

Parameter	Value
Product Name	TERT Rabbit pAb
Remarks/Alias	TP2, TRT, CMM9, EST2, TCS1, hTERT, DKCA2, DKCB4, hEST2, PFBMFT1, TERT
Species	Human
GeneID (Human)	7015
GeneID	7015
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1050-1132 of human TERT (NP_937983.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O14746
Protein Weight	127kDa
Shipping	Ice bag