

Anti-Tes Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Ugi8

Item Number: CM0028480



Introduction

Rabbit polyclonal antibody against human TES (Testin), validated for WB, IHC-P, ELISA. Reacts with human and rat. TES is a scaffold protein involved in cell adhesion and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	TES Rabbit pAb
Remarks/Alias	TESS; TESS-2; TES
Species	Human
Gene ID (Human)	26136
Gene ID	26136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TES (NP_056456.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UGI8
Protein Weight	48kDa
Shipping	Ice bag