

Anti-Tet1 Monoclonal Antibody For Wb And Elisa - Q8Nfu7

Item Number: CM0007695



Introduction

High-quality rabbit monoclonal antibody targeting human TET1, suitable for Western blot and ELISA applications. Reacts with human, mouse, and rat samples, ideal for epigenetic and developmental research.

[Learn More](#)

Parameter	Value
Product Name	TET1 Rabbit mAb
Remarks/Alias	LCX; CXXC6; bA119F7.1; TET1
Species	Human
GeneID (Human)	80312
GeneID	80312
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFU7
Protein Weight	235kDa
Shipping	Ice bag