

Anti-Tet1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q8Nfu7

Item Number: CM0018478



Introduction

Validate TET1, a key DNA demethylation and chromatin regulator, with this rabbit polyclonal antibody. Tests in WB, IF/ICC, IHC-P, ELISA; reactive across human, mouse, rat samples. (Q8NFU7)

[Learn More](#)

Parameter	Value
Product Name	TET1 Rabbit pAb
Remarks/Alias	LCX; CXXC6; bA119F7.1; TET1
Species	Human
Gene ID (Human)	80312
Gene ID	80312
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFU7
Protein Weight	235 kDa
Shipping	Ice bag