

Anti-Tet1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q3UrK3

Item Number: CM0024563

Introduction

Rabbit polyclonal antibody against mouse TET1, validated for WB, IHC-P, and ELISA. Detects a ~219 kDa band; cross-reacts with human, mouse, and rat. Suitable for DNA demethylation and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	TET1 Rabbit pAb
Remarks/Alias	LCX; Cxxc6; mKIAA1676; D10ErtD17e; 2510010B09Rik; TET1
Species	Mouse
GeneID	52463
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1749-1902 of mouse TET1 (NP_001240786.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3URK3
Protein Weight	219kDa
Shipping	Ice bag