

Anti-Tet2 Rabbit Pab For Wb, Ihc-P, Elisa - Q6N021

Item Number: CM0018006



Introduction

TET2 Rabbit pAb, a polyclonal antibody validated for WB, IHC-P, and ELISA, recognizes human TET2 (Q6N021) and cross-reacts with mouse and rat. Suitable for DNA methylation and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	TET2 Rabbit pAb
Remarks/Alias	MDS; IMD75; KIAA1546; TET2
Species	Human
Gene ID	54790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1833-2002 of human TET2 (NP_001120680.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6N021
Protein Weight	224kDa
Shipping	Ice bag