

Anti-Tet2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6N021

Item Number: CM0026502



Introduction

Validated for WB, IHC-P, IF/ICC, and ELISA, this rabbit polyclonal antibody targets human TET2 (Q6N021), a methylcytosine dioxygenase involved in DNA demethylation, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	TET2 Rabbit pAb
Remarks/Alias	MDS; IMD75; KIAA1546; TET2
Species	Human
GeneID (Human)	54790
GeneID	54790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1833-2002 of human TET2 (NP_001120680.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6N021
Protein Weight	224kDa
Shipping	Ice bag