

Anti-TEX101 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9By14

Item Number: CM0022682



Introduction

Rabbit polyclonal antibody targeting human TEX101. Validated for WB, IHC-P, ELISA; reacts with human, mouse, and rat samples. Ideal for fertilization and reproductive biology studies.

[Learn More](#)

Parameter	Value
Product Name	TEX101 Rabbit pAb
Remarks/Alias	SGRG; CT131; GTPR867; NYD-SP8; PRO1884; SPATA44; TES101RP; TEX101
Species	Human
GeneID	83639
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-150 of human TEX101 (NP_113639.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY14
Protein Weight	27kDa
Shipping	Ice bag