

Anti-Tfam Rabbit Monoclonal Antibody For Wb - Q00059

Item Number: CM0007434



Introduction

Rabbit monoclonal antibody against human TFAM, validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Ideal for mitochondrial transcription and mtDNA maintenance research.

[Learn More](#)

Parameter	Value
Product Name	TFAM Rabbit mAb
Remarks/Alias	TCF6; MTTF1; MTTFA; TCF6L1; TCF6L2; TCF6L3; MTDP515; TFAM
Species	Human
GeneID (Human)	7019
GeneID	7019
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human
SWISS	Q00059
Protein Weight	29 kDa
Shipping	Ice bag