

Anti-Tfam Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q00059

Item Number: CM0020160



Introduction

High-quality rabbit polyclonal antibody against human TFAM, suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat, ideal for mitochondrial transcription and mtDNA maintenance research.

[Learn More](#)

Parameter	Value
Product Name	TFAM Rabbit pAb
Remarks/Alias	TCF6; MTTF1; MTTFA; TCF6L1; TCF6L2; TCF6L3; MTDP515; TFAM
Species	Human
GeneID (human)	7019
GeneID	7019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TFAM (NP_003192.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00059
Protein Weight	29kDa
Shipping	Ice bag