

## Anti-Tfam Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q00059

Item Number: CM0018266



### Introduction

Rabbit polyclonal antibody targeting human TFAM, suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial transcription and mtDNA maintenance research.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | TFAM Rabbit pAb                                                                                      |
| Remarks/Alias    | TCF6; MTTF1; MTTFA; TCF6L1; TCF6L2; TCF6L3; MTDPS15; TFAM                                            |
| Species          | Human                                                                                                |
| Gene ID (Human)  | 7019                                                                                                 |
| Gene ID          | 7019                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TFAM (NP_003192.1) |
| Source           | Rabbit                                                                                               |
| Category         | Polyclonal Antibodies                                                                                |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                    |
| SWISS            | Q00059                                                                                               |
| Protein Weight   | 29 kDa                                                                                               |
| Shipping         | Ice bag                                                                                              |