

Anti-Tfap2B Polyclonal Antibody For Wb, Elisa - Q92481

Item Number: CM0029969



Introduction

Rabbit polyclonal antibody targeting human TFAP2B, a transcription factor involved in development. Suitable for Western blot and ELISA. Specific to human. Predicted molecular weight 50 kDa.

[Learn More](#)

Parameter	Value
Product Name	TFAP2B Rabbit pAb
Remarks/Alias	PDA2; AP-2B; AP2-B; AP-2beta; TFAP2B
Species	Human
GeneID (human)	7021
GeneID	7021
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TFAP2B (Q92481).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92481
Protein Weight	50kDa
Shipping	Ice bag