

Anti-Tfap2C Polyclonal Antibody For Wb, Elisa - Q92754

Item Number: CM0025928

Introduction

Rabbit polyclonal antibody against human TFAP2C (Transcription factor AP-2 gamma), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for developmental biology, cancer research, and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	TFAP2C Rabbit pAb
Remarks/Alias	ERF1; TFAP2G; hAP-2g; AP2-GAMMA; TFAP2C
Species	Human
GeneID (Human)	7022
GeneID	7022
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of human TFAP2C (NP_003213.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92754
Protein Weight	49kDa
Shipping	Ice bag