


Anti-Tfb1M Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wvm0

Item Number: CM0015812



Introduction

Rabbit polyclonal antibody targeting human TFB1M (dimethyladenosine transferase 1, mitochondrial) for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying mitochondrial ribosome assembly and mtDNA transcription.

[Learn More](#)

Parameter	Value
Product Name	TFB1M Rabbit pAb
Remarks/Alias	CGI75; mtTFB; CGI-75; mtTFB1; TFB1M
Species	Human
Gene ID (Human)	51106
Gene ID	51106
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-346 of human TFB1M (NP_057104.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVM0
Protein Weight	40kDa
Shipping	Ice bag