

Anti-Tfcp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q12800

Item Number: CM0030792

Introduction

Rabbit polyclonal antibody against alpha-globin transcription factor CP2 (TFCP2/LSF/SEF). Validated for WB and ELISA. Reacts with human and mouse. Suitable for chromatin remodeling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	TFCP2 Rabbit pAb
Remarks/Alias	LSF; SEF; LBP1C; LSF1D; TFCP2C; TFCP2
Species	Human
GeneID (Human)	7024
GeneID	7024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TFCP2 (NP_005644.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12800
Protein Weight	57kDa
Shipping	Ice bag